

estheticpro

The EstheticPro Micropen is versatile, powerful, adjustable, efficient, and economical: making it an ideal tool for an aesthetician. It can be used for microneedling and nanoinfusion treatments.

Features:

- ☀ 5 Adjustable Speeds - 8,000 rpm's to 18,000 rpm's
- ☀ Battery Powered
- ☀ 0 to 2.0mm adjustable needle depth
- ☀ 4 Different Cartridge Choices: (single use only)
 - Silicone Nano* - Safe for all skin types including sensitive skin or metal allergies
 - Surgical Stainless Steel Nano* - Most popular
 - 36 Needle* - Aesthetic grade cartridge
 - 12 Needle* - Medical grade cartridge

Power: Lithium ion battery is charged via DC adapter to AC outlet. Best to charge the pen in smaller bursts (i.e. 2 hours) than leaving it plugged in for long periods of time (i.e. over night). To extend the battery life, charge the pen at least once a month.

Cleaning: Wipe with disinfectant before and after each use. Do not submerge in liquid disinfectant.

2 Year warranty

